

Data sheet

ultraSTRETCHER "0710"

Stretcher according to DIN 13024 K made of ultraCHEM fabric.

ultraCHEM: Chemical fiber fabric, Highly tear-resistant and impactresistant, Waterproof disinfectable,

2x foldable, with 4 slide feet.

Robust, easy, proven

All stretchers own:

- Small stowage / pack size
- Aluminum base frame
- Integrated pillow bag
- Solid lockable steel hinges for folding
- Two adjustable straps with quick closure for secure fixation of the injured.

SKU: SAN-0104-4G

Weight: 7.5 kg

Colour: Gray

EAN: 4030281000792

DIMENSIONS

Width: 556 mm

Height: 137 mm

Length: 2302 mm

