



RESCUE STRETCHERs



PERFECT SOLUTIONS



FOR THE MOST DIVERSE FIELDS
OF APPLICATION



The specific rescue situation and the constitution of the patient determine the choice of suitable rescue equipment. The correct choice of rescue stretcher is decisive for a successful patient transport.

Whether for water rescue, X-ray capable or for the rescue of obese patients:

the ultraMEDIC rescue stretcher series offers the perfect solution for different fields of application.



ultraSCOOP STRETCHER & ultraSCOOP STRETCHER X-RAY

SCOOP STRETCHER

SCOOP STRETCHER



SAN-0125



SAN-0288



6,5 kg



1150 x 430 x 90 mm smallest size
1650-2010 x 430 x 90 mm max. extended



incl. 3 patient safety belts



ultraSCOOP STRETCHER
SAN-0125



7,2 kg



1180 x 450 x 80 mm smallest size
1650-2010 x 450 x 80 mm max. extended



X-ray-compatible



incl. 3 patient safety belts



ultraSCOOP STRETCHER X-RAY
SAN-0288

MATERIAL

Colorless anodized special aluminum alloy

AREA OF APPLICATION

For gentle and safe lifting and repositioning of severely injured persons. For this purpose, the scoop halves are opened and pushed sideways under the patient.

RECOMMENDED

In case of suspected spinal injuries

MATERIAL

Plastic on aluminum frame

The plastic surface is closed and prevents the penetration of liquids.

AREA OF APPLICATION

For gentle and safe lifting and repositioning of severely the scoop halves are opened and pushed sideways under the patient.

RECOMMENDED

In case of suspected spinal injuries

The entire spinal area of the scoop stretcher is **radiolucent**.



info@ultramedic.de



ultraMEDIC GmbH
Im Bruch 11-15
56567 Neuwied



www.ultramedic.de

ultraSCOOP STRETCHER CARBON & ultraCOMBI-STRETCHER II

SCOOP STRETCHER

COMBI-DEVICE



SAN-0297



SAN-0298



5,6 kg



1180 x 450 x 80 mm smallest size
1650-2010 x 450 x 80 mm max. extended



X-ray-compatible



incl. 3 patient safety belts



ultraSCOOP STRETCHER CARBON
SAN-0297



8 kg



1870 x 420 x 57 mm



Suitable for X-rays and CT



incl. 3 patient safety belts



ultraCOMBI-STRETCHER II
SAN-0298

MATERIAL

Ultra light thanks to carbon fibers

AREA OF APPLICATION

For gentle and safe lifting and repositioning of severely injured persons. For this purpose, the bucket parts are opened and pushed sideways under the patient.

RECOMMENDED

In case of suspected spinal injuries

MATERIAL

Lightweight, durable PEHD outer material with carbon core.

AREA OF APPLICATION

Due to the double-sided shackle locks and the particularly flat edge, the shoveling and locking of the shovel halves is quick and easy. The concave shape allows completely immobilization for easy and safe transport. The gentle treatment of the patient without subsequently lifting the head or legs distinguishes this model from other combination devices.

RECOMMENDED

Combines all the advantages of a scoop stretcher and a spine-board and due to its high loading capacity of max. 300 kg and flat design, it is ideally suited for the transport of obese patients.



info@ultramedic.de



ultraMEDIC GmbH
Im Bruch 11-15
56567 Neuwied



www.ultramedic.de

ultraCOMBI-STRETCHER & S-BOON

COMBI-DEVICE

SPINEBOARD



SAN-0296

- 8,1 kg
- 1860 x 430 x 70 mm
- Suitable for X-rays and CT
- incl. 3 patient safety belts

Art. No. ultraCOMBI-STRETCHER
SAN-0296

MATERIAL

Plastic surface is insensitive, easy to clean and disinfect.

AREA OF APPLICATION

The ultraCOMBI-STRETCHER combines in its characteristics all the advantages of a scoop stretcher and a spineboard. Thus it offers the optimal combination for rescuing, transporting and repositioning the patient.

RECOMMENDED

Universally applicable



SAN-0294

- 7,5 kg
- 1850 x 430 x 50 mm
- MRI-compatible
- Suitable for X-rays and CT
- Buoyant (buoyancy approx. 20 kg / 21 kN)
- SAN-0295 incl. 3 patient safety belts



Art. No. S-BOON
SAN-0295 with straps
SAN-0294 without straps

MATERIAL

A frame construction originating from aircraft construction with integrated **carbon fiber inlay** or **carbon core**.

MORE THAN A BOARD

Slightly concave construction, allowing sufficient space for even "reaching under" and lifting is possible -> **use possible with working gloves**.

AREA OF APPLICATION

It is mainly used for lifting and repositioning patients with severe head and spinal injuries.
Extremely resilient.

RECOMMENDED

Ideally suited for water rescue operations



info@ultramedic.de



ultraMEDIC GmbH
Im Bruch 11-15
56567 Neuwied



www.ultramedic.de

our
Responsibility

RESCUE STRETCHER_s

There is undoubtedly one thing in life where there can be no compromise: when it comes to helping people in emergency and borderline situations! For this, ultraMEDIC is the partner at your side. Our name stands for highly specialized quality products for **FIRST AID, EMERGENCY and RESCUE.**

We develop and produce our own solutions at our headquarter in Neuwied / Rhein with passion and many years of know-how. Know-how for all those who help, rescue and recover: fire departments, rescue services, occupational health and safety. ultraMEDIC is probably the only manufacturer of stretchers in Germany with the predicate "Made in Germany". Technical development, production and support are united under one roof. The intensive exchange of experience with users and the rapid and flexible implementation of individual requirements and wishes in innovative samples and prototypes. Prototypes are a matter of course for us. Because we never lose sight of one thing: Your mission is our responsibility. - Your demand is our challenge! You can rely on solutions from ultraMEDIC in any situation and focus on the essentials:

REDUCE DISTRESS, ALLEVIATE PAIN, SAVE PEOPLE!



info@ultramedic.de



ultraMEDIC GmbH
Im Bruch 11-15
56567 Neuwied



www.ultramedic.de

ACCESSORIES

ultraHEAD-LOCK III Head fixation set



Art. No. SAN-0291

kg 0,52 kg

mm 350 x 230 x 140 mm

Icon Suitable for MRI, X-rays and CT

Icon Buoyant

+ COMPATIBLE with
SAN-0298
ultraCOMBI-STRETCHER II

ultraHEAD-LOCK II Head fixation set



Art. No. SAN-0285
SAN-0285-P Replacement pads

kg 0,6 kg

mm 380 x 220 x 170 mm

Icon incl. storage bag SAN-0160

Icon Buoyant

Icon Suitable for MRI, X-rays and CT

+ COMPATIBLE with
SAN-0294 / SAN-0295
S-BOON Spineboard
SAN-0296
ultraCOMBI-STRETCHER

ultraHEAD-LOCK I Head fixation set



Art. No. SAN-0284

kg 1,0 kg

mm 430 x 280 x 160 mm

Icon Suitable for MRI, X-rays and CT

Icon Buoyant

+ COMPATIBLE with
SAN-0125
ultraSCOOP STRETCHER
SAN-0288
ultraSCOOP STRETCHER X-RAY
SAN-0296
ultraCOMBI-STRETCHER
SAN-0297
ultraSCOOP STRETCHER CARBON
SAN-0294 / SAN-0295
S-BOON Spineboard

ultraBELT SET Replacement straps set



Art. No. SAN-0084-K

kg 0,9 kg

Icon 3 Replacement straps incl. storage bag

+ COMPATIBLE with
SAN-0125
ultraSCOOP STRETCHER
SAN-0288
ultraSCOOP STRETCHER X-RAY
SAN-0297
ultraSCOOP STRETCHER CARBON
SAN-0298
ultraCOMBI-STRETCHER II
SAN-0294 / SAN-0295
S-BOON Spineboard

ultraSTRAPS Restraint belt system with Velcro



Art. No. SAN-0281

kg 0,4 kg

Icon Suitable for MRI, X-rays and CT

ultraSTRAPS Xtreme Restraint belt system with clip lock



Art. No. SAN-0286

kg 1 kg

Icon Suitable for MRI, X-rays and CT

+ COMPATIBLE with
SAN-0294
S-BOON Spineboard