



NO 2nd CHANCE for FIRST AID



EMERGENCY & FIRST AID

IN THE COMPANY

DEFIBRILLATOR HEARTSINE PAD 360P



Fully automatic defibrillator for use with adults and children. The device analyzes the heart rhythm and, if necessary, automatically delivers an electric shock to save the patient.

SAN-8620

200 x 184 x 48 mm

DEFIBRILLATOR LIFEPAK CR2



Fully automatic, intuitive defibrillator for use with adults and children. A daily self-test gives you the assurance that the device is always ready for use. The integration into a WLAN also facilitates diagnosis for the emergency response team.

SAN-8581

226 x 274 x 97 mm

FIRST AID BACKPACK ultraRESCUE LITE



Light and compact rescue backpack with extended first aid filling for use in the factory or on the construction site. Work like the pros - 4 colored module pockets in the main compartment provide for order and a fast access to the contents. In the front pocket there is even room for a defibrillator.

SAN-400-R-FA

450 x 270 x 180 mm



info@ultramedic.de



ultraMEDIC GmbH
Im Bruch 11-15
56567 Neuwied



www.ultramedic.de



No 2ND CHANCE FOR FIRST AID

EMERGENCY & FIRST AID IN THE COMPANY

FIRST AID EQUIPMENT IN COMPANIES

First Aid Kit (DGUV Rule 100-001 - May 2014 edition)

Administrative and commercial operations

1-50



51-300



301-600



+ per 300+



Manufacturing, processing and comparable companies

1-20



21-100



101-200



+ per 100+



Construction sites and construction site-like operations

1-10



10-50



51-100



+ per 50+



Couches (Recommendation)

up to 20 1x

up to 300 4x

up to 1000

up to 50 2x

up to 500 5x

7x

up to 100 3x

up to 750 6x

Industry sectors First Aid Kit



First aid kit with a wall bracket, large, basic contents in accordance with DIN 13157 + industry-specific First Aid products. Industry-specific cases BAU and WERK also contain a medical replant bag, among other things.

- Industry-specific pictogram
- Also available for other industries

445 x 340 x 150 mm

- SAN-0175-BAU
- SAN-0175-WERK
- SAN-0175-OFFICE

First Aid Box BASIC



First aid kit with a wall bracket, small, with filling according to DIN 13157.

- SAN-0170-OR
- 285 x 215 x 120 mm

First Aid Box SUPER II



First aid kit with a wall bracket, medium-sized, with filling according to DIN 13169.

- SAN-0171-OR-G
- 330 x 250 x 130 mm



No 2ND CHANCE FOR FIRST AID

EMERGENCY & FIRST AID IN THE COMPANY

ultraTRAFFIC vehicle first aid box



Vehicle first aid box
made of ABS plastic



SAN-0176
250 x 170 x 80 mm

ultraTRAFFICBAG



Vehicle first aid nylon bag



SAN-0176-T
250 x 170 x 80 mm

ultraKIT



First Aid Kit with zipper

SAN-0176-K1
content for travelling 170 x 125 x 70 mm

SAN-0176-K2
with content 110 x 80 x 40 mm

ultraWORKS



Plastic first aid box,
filled according to
DIN 13157.



SAN-0165
250 x 170 x 80 mm

Standing cabinet ultraLOCKER 1350



Narrow cabinet made of quality thin sheet metal,
with all-round door seal and an emergency key
compartment with glass cover, can be ordered
with or without filling assortment. Filled cabinets
contain:

- Stretcher
- Blanket
- Rescue seat
- First aid kit (DIN 13157)
- Rescue and evacuation blanket
- Flashlight with red/green color filter

300 x 200 x 2000 mm

SAN-0197-OR (without content)

SAN-0197-1OR (with content)

First Aid cabinet ultraTOP-LOCKER 023



Add-on first-aid cabinet with filling
range according to DIN 13169,
with 2 intermediate shelves

SAN-0069-OR
490 x 200 x 560 mm

Rest couch ultraREST 200



2-part rest couches with adjustable head section,
extra thick foam pad (ÖkoTex Standard 100),
fabric-reinforced artificial leather, easy to clean, lying
surface height approx. 500 mm, available in 15
colors.

- | | | |
|---------------|--------------|------------|
| Black | Melon | Dark Blue |
| White | Medium Brown | Fir Green |
| Light Beige | Red | Petrol |
| Sahara | Burgundy | Light Gray |
| Pastel Yellow | Medium Blue | Dark Gray |

SAN-0147-DBL
1900 x 700 x 500 mm



No 2ND CHANCE FOR FIRST AID

EMERGENCY & FIRST AID IN THE COMPANY

ultraTHERM protective cover



Optimal thermal insulation in case of emergency.
Can be used as a blanket or complete cover: Velcro fasteners & side zippers allow flexible application.
Washable & disinfectable.

 SAN-7072

 750 x 2200 mm (unfolded)

ultraTHERM patient blanket



Bacteriostatic blanket for warming or keeping warm patients.
Washable & disinfectable.
Polyester fiber filling.

 SAN-7071


 1100 x 1900 mm

ultraRESCUE-CHAIR



Rescue and evacuation chair for transporting patients in emergencies. Comfortable rolling on flat surfaces.
Fold-out handles allow carrying by 2 persons in uneven terrain or staircases.

 SAN-9201

 510 X 920 X 600 mm

ultraRESCUE-CHAIR STAIRS



Rescue and evacuation chair for transporting patients in emergencies. Comfortable rolling on flat surfaces. Fold-out handles and extendable transport handle allow the chair to be carried by 2 persons on uneven terrain or staircases.
Braking system for sliding down stairs.

 SAN-9202

 540 X 1000 X 720 mm



info@ultramedic.de



ultraMEDIC GmbH
Im Bruch 11-15
56567 Neuwied



www.ultramedic.de





No 2ND CHANCE FOR FIRST AID

EMERGENCY & FIRST AID IN THE COMPANY

Basket Stretcher with lifting bridles





Robust rescue basket for rescue and transport of injured persons, e.g. from scaffolds.
The stretcher can be used horizontally and vertically and due to the stable construction, it can also be used as a descender cage.

-  SAN-0087-1 (Basket stretcher)
-  SAN-0086 (Lifting bridles)

Patient stretcher ultraSTRETCHER 0350



Stretcher with 4 sliding feet according to DIN 13024 N.
1x longitudinally foldable, made of highly tear-resistant man-made fiber fabric, resistant, waterproof, disinfectable.

-  SAN-0101-4G
-  2310 x 550 x 140 mm (ready for use)
1910 x 140 x 140 mm (folded)



info@ultramedic.de



ultraMEDIC GmbH
Im Bruch 11-15
56567 Neuwied



www.ultramedic.de



No 2ND CHANCE FOR FIRST AID

EMERGENCY & FIRST AID IN THE COMPANY

Disinfectant Dispenser



- ⊕ Contactless removal
- ⊕ Wall mounting or Mounting on free-standing column
- ⊕ Refill container with 1l volume
- ⊕ Battery or mains operation possible (*power cord can be ordered separately*)
- ⊕ 2 dispensing quantities adjustable:
 - 3ml (*recommended quantity for hand disinfection*)
 - 2x 3ml



- SAN-0050-2 ultraTOWER II
- SAN-0050-2-S ultraDISPENSER

BREASY MNS IIR Filter Mask



Certified surgical mask from the group's own production at the Neuwied site.

CE Made in Germany

□ SAN-0028-IIR

BREASY FFP2 Folding Mask



Certified FFP2 mask from the group's own production at the Neuwied site.

CE ₂₁₆₃ Made in Germany

□ SAN-0028-FFP2

ultraMASK



Respiratory mask with oxygen connection, integrated, exchangeable bacteria filter and one-way valve, 2 alcohol pads, 2 disposable gloves. Delivery in EVA pocket box with practical carrying loop.

■ SAN-0256-P

■ 140 × 113 × 50 mm

ultraMEDIC

There is undoubtedly one thing in life where there can be no compromise: when it comes to helping people in emergency and border situations! For this, ultraMEDIC is the partner at your side. Our name stands for highly specialized quality products for **First Aid, Emergency and Rescue**.

We develop and produce our own solutions at our headquarter in Neuwied/Rhein with passion and many years of know-how. know-how for all those who help, rescue and recover: fire departments, rescue services, occupational health and safety. ultraMEDIC is probably the only manufacturer of stretchers in Germany with the predicate "Made in Germany". Technical development, production and support are united under one roof. The intensive exchange of experience with users and the rapid and flexible implementation of individual requirements and wishes in innovative samples and prototypes. Prototypes are a matter of course for us. Because we never lose sight of one thing: Your mission is our responsibility.- Your demand is our challenge!

You can rely on solutions from ultraMEDIC in any situation and focus on the essentials:

REDUCE DISTRESS, ALLEVIATE PAIN, SAVE PEOPLE!