



RED has a SIGNAL EFFECT

HIGHLIGHTS

FROM THE
BACKPACK & BAG RANGE



RESCUE I
RESCUE BACKPACK

A TRIED-AND-TESTED CLASSIC WITH
PRICE-PERFORMANCE GUARANTEE



RESCUE I PLUS
RESCUE BACKPACK WITH
1 ADDITIONAL BAG

THE PLUS OF EXTRA SPACE FOR
ACCESSORIES, AED OR PERSONAL
EQUIPMENT



ultraRESCUE LITE
RESCUE BACKPACK
COMPACT

LITTLE SPACE, LOTS OF MATERIAL?
YOUR SOLUTION - OUR COMPACT
SPACE-SAVING BACKPACK.



ultraDOC ROT
DOCUMENT POUCH

MORE THAN JUST A
DOCUMENT POUCH -
THE DIFFERENCE IS IN THE
DETAILS.



RESCUE I & RESCUE I PLUS Rescue backpacks

PRICE meets PERFORMANCE



The tried-and-tested and versatile **RESCUE I** now has a big brother: the **RESCUE I PLUS**. Just as stable and robust, and now with an extra front pouch for storing extra equipment.

RESCUE I
Part-No. **SAN-500-R**

RESCUE I PLUS
Part-No. **SAN-500-R-Z**

MULTI-TALENT



RESCUE I



FIRST
RESPONDER
DIN 13155



DOOR
OPENING



FIRE &
RESCUE
DIN 14142



RESCUE I PLUS



On the inside, the RESCUE I stays flexible yet tidy thanks to colour-coded modular compartments incl. a holder for an oxygen cylinder.



info@ultramedic.de



ultraMEDIC GmbH
Im Bruch 11-15
56567 Neuwied



www.ultramedic.de

ultraRESCUE LITE Compact rescue backpack



WORKING like a PRO

With external dimensions of 45 x 27 x 21 cm, this compact rescue backpack is ideal for a wide variety of applications and leaves nothing to be desired when it comes to storage. 4 outer pouches and 4 colour-coded modular compartments inside can hold a lot of material. And if necessary, the main compartment can be extended by a further 5 cm - with additional space for a cervical support.

Part-No. **SAN-400-R**

STRUCTURE meets TIDINESS



Colour coded

Extension --
Red zip



Organising systems

in all
compartments



Identification pocket window

Opening downwards
to protect from
splash water



Handle area

Reinforced
handle and
material



Professional

Modular bags from
the professional
backpack
ultraRESCUE I



VARIANTS

SAN-400-R-L

Backpack without
modular bags

SAN-400-R

Backpack with
4 modular bags



SAN-400-R-FR

Backpack incl. DIN 14142
+ special equipment
„helpers on site“



SAN-400-R-FA

Backpack incl. DIN 13157
+ special equipment
„Company“



SAN-400-R-SP

Backpack incl. DIN 13157
+ special equipment
„Sport“



info@ultramedic.de



ultraMEDIC GmbH
Im Bruch 11-15
56567 Neuwied



www.ultramedic.de

ultraDOC ROT Document pouch



ORGANIZED in ACTION

The robust and washable document pouch is equally versatile and practical. Apart from a large number of pockets for documents, cards etc., a clipboard and a tablet holder are also integrated. Using velcro, the clipboard can also be secured to the outside of the bag for use as a writing pad.

Part-No. **SAN-7225-R**

THE MINI OFFICE



Everything has its place inside: documents, pens, instructions, mission cards, personal equipment and even individual keys.



The clipboard can be secured to the outside using a velcro system and when the zip is closed, the ultraDOC serves as a writing pad which can be comfortably strapped around the shoulders.



On the back, documents can be stored visibly and protected from the rain or mission scenarios can be outlined on insert cards.



info@ultramedic.de



ultraMEDIC GmbH
Im Bruch 11-15
56567 Neuwied



www.ultramedic.de